

As we recognize the important role of early intervention in cognitive development, more and more children use mobility bases at an earlier age. VARILITE JUNIOR seating and positioning products were designed with the knowledge that kids are not small adults—that's why we focused on the not-so-small details that make our products just for them.

The VARILITE JUNIOR line features a full-range of size-appropriate seating and positioning items to support children in their mobility device, allowing them to go where their imagination takes them. By exploring their world, children can gain confidence, increased peer interaction and independence. And VARILITE JUNIOR seating systems are there to support them along the way.



#### Sized For Children

VARILITE JUNIOR products were designed with children in mind using anthropometric data.



#### Lightweight

VARILITE JUNIOR products combine lightweight positioning materials to support kids on the move.



#### Designed for Growth

VARILITE JUNIOR products allow for growth with simple-to-adjust hardware and easy-to-fit soft goods.



#### Low Maintenance for Busy Parents

VARILITE JUNIOR products are simple to use to keep up with busy families.

VARILITE JUNIOR Seating System

## VARILITE JUNIOR Secondary Supports

VARILITE Secondary Supports are the perfect accessory to our JUNIOR Seating System. Lightweight yet durable, our chest harnesses, hip belts and Band-It™ Ankle Supports provide gentle support and positioning. Machine washable.

### JUNIOR Contoured Chest Harness

VARILITE Contoured Chest Harness is designed to achieve an optimal fit using soft molded pads and our Dynamic Positioning Strap.

Part #	Description	Pad Size	Webbing
87115	Side squeeze, shoulder cam, front pull*	9x1.5 in. (22x3.8 cm)	3/4 in. (1.9 cm)

\*Bottom Pull Harness Accessory available



### JUNIOR Hip Belts

VARILITE Hip Belts enhance posture and function through anterior pelvic support. Available in a wide range of styles to meet a variety of positioning needs.

Part #	Description	Buckle Material	Pad Size	Webbing
83012	2-point, side squeeze, single pull	metal	4x1.5 in. (10x3.8 cm)	1 in. (2.54 cm)
85012	2-point, side squeeze, single pull	plastic	4x1.5 in. (10x3.8 cm)	1 in. (2.54 cm)
83016	2-point, side squeeze	metal	4x1.5 in. (10x3.8 cm)	1 in. (2.54 cm)
85016	2-point, side squeeze	plastic	4x1.5 in. (10x3.8 cm)	1 in. (2.54 cm)
85014	4-point, side squeeze	metal	4x1.5 in. (10x3.8 cm)	1 in. (2.54 cm)



### JUNIOR Band-It™ Ankle Supports

VARILITE Band-It™ Ankle Supports position lower extremities which in turn help position the entire body. Our Dynamic Positioning Strap allows freedom of ankle movement.

Part #	Description	Ankle Circumference
89003	Side squeeze	4-6 in. (10-15 cm)



Pediatric Seating and Positioning System



A Better Way to Sit™



A Better Way to Sit™



## VARILITE JUNIOR Seat Cushion

The JUNIOR Seat Cushion provides a comfortable, stable foundation for exploration and adventure. Moderate skin protection and positioning is achieved using VARILITE air-foam floatation technology and a soft tailored cover.

### Features



**Multi-stiffness foam** and a non-intrusive medial thigh separator promotes immersion, stability and positioning.



**Water resistant, washable covers** and a waterproof inner cushion let kids be kids.



**Quick adjustments** can be made using the self-inflating valve.

Part #	Description	Height	Weight	Dimensions
77000110102	JUNIOR Seat Cushion	2.9 in. (7.4 cm)	9.2 oz. (261 g)	10x10 in. (25x25 cm)

## VARILITE JUNIOR Back Support

Lightweight and easy to use, the JUNIOR Back Support is perfect for children and families on the go. Moderate positioning is achieved using an adjustable VARILITE air-foam floatation cushion and versatile mounting hardware. Offered in mid and deep contoured styles.


### Features



**Compact, minimal design** directs focus to the child, not the equipment.



**Simple adjustments** allow the back system to grow with your child.




**Lightweight** air-foam construction coupled with lightweight aluminum allows your child to manually propel with ease.



**Simple operation** means less time fussing with hardware and more time on the go.

## JUNIOR Back Support Technical Information



	Mid Back	Deep Back
A. Back height*	9.25 in. (23.5 cm)	9.25 in. (23.5 cm)
B. Bottom of back to beginning of vertical slot	1.0 in. (2.5 cm)	1.0 in. (2.5 cm)
C. Bottom of back to top of vertical slot	7.0 in. (17.8 cm)	7.0 in. (17.8 cm)
D. Range of hardware adjustment in vertical slot	4.91 in. (12.5 cm)	4.91 in. (12.5 cm)
E. Outside width of shell	8.6 in. (21.8 cm)	9.9 in. (25.1 cm)
F. Inside width of support surface	7.6 in. (19 cm)	8.9 in. (22.6 cm)
G. Height of fixed lateral support	N/A	5.7 in. (14.5 cm)
H. Height of lateral slot	3 in. (7.6 cm)	N/A
Fits wheelchairs sizes	9-11 in. (23-28 cm)	9-11 in. (23-28 cm)
Weight, including hardware	2.4 lbs. (1.1 kg)	2.9 lbs. (1.3 kg)
Part #	4430111020	4430411020

\*Shell only. Complete back system height may vary by 1 in. (2.5 cm) due to the cushion rollover.