

# The ultimate mobility tool

**Batec handbikes** are a comprehensive solution to your mobility problems:

- ✓ Go from reduced to **augmented mobility**; from needing help to being **independent** and **proactive**.
- ✓ Get a **2 in 1**, without giving up your manual wheelchair.
- ✓ **Overcome obstacles** and **drive over bumpy terrains**.
- ✓ Attach our handbikes **in just seconds** and **with no need for transfers**.
- ✓ Also available for **quadriplegics** and **hemiplegics**.



You CAN get used  
to being in a  
**wheelchair**



But not to giving up on

**LIFE**

[www.batec-mobility.com](http://www.batec-mobility.com)



Distributed in Canada by  
49 Bespoke Inc.  
51 Caldari Road Units 6-7B,  
Concord, Ontario L4K 4G3  
[www.49bespoke.com](http://www.49bespoke.com)  
Ph: +1.866.827.4548  
Fax: +1.877.731.7745  
Email: [49@49bespoke.com](mailto:49@49bespoke.com)



**BATEC ELECTRIC 2**  
The perfect balance



# BATEC ELECTRIC 2, the perfect balance

The **BATEC ELECTRIC** is the latest in our saga of handbikes. It features the latest **Batec Mobility** technology and an **18" Wheel** designed exclusively for this model, which means you can use it both in **urban** and **off-road settings**.



All-in-one LCD screen with gear selector.



18" Wheel, ultra-grip tyre and aluminium double-wall rim.



250 RPM 900W high torque brushless motor with reverse gear. Maximum speed 25 Km/h.



Avid BB5 mechanical brakes, aluminium levers and parking brake.



700-lumen front light with position lights / high beam options. Dual LED tail light on the stand frame.



Horn and USB port.



Removable right dropout for easy inner tube or tyre change.



Removable QR weights.



Autonomy of 22 km \* with a 36v - 280 Wh battery. Autonomy of 40 km \* with a 36v - 522 Wh battery. 4A Charger.




Matte black aluminium handlebar with folding system for convenience.



G&B handlebar option without levers for quadriplegics.



Handlebar option with right or left Monolever brakes for hemiplegics.



## COLOURS



\* Autonomy calculated under the following conditions: 75 - kg user, ultra-light wheelchair, flat terrain, no wind and steady speed of 15 km/h. May vary +/- 20% depending on the battery. Number/version: MM-008/3.0